


Excelsior Laser Head Dimensions

Model	532 nm (Multi Mode Non-collimated)	375 nm, 405 nm, 440 nm, 473 nm, 488 nm, 532 nm (Single Mode, Multi Mode Collimated), 542 nm, 561 nm, 642 nm, 785 nm, 1064 nm	Direct Diode and DPSS CDRH models	505 nm, 515 nm, 594 nm, 561 nm (100 & 150 mW), 532 nm (300 mW), DPSS XC CDRH models
A Length	3.33 in (84.5 mm)	3.74 in (95 mm)	4.33 in (110 mm)	4.53 in (115 mm)
B Width	1.1 in (28 mm)	1.1 in (28 mm)	1.65 in (42 mm)	1.97 in (50 mm)
C Height	1.44 in (36.5 mm)	1.44 in (36.5 mm)	1.44 in (36.5 mm)	1.71 in (43.5 mm)
D Beam Height	0.75 in (19 mm)	0.75 in (19 mm)	0.75 in (19 mm)	0.75 in (19 mm)
E Beam Exit (from side)	0.55 in (14 mm)	0.55 in (14 mm)	0.83 in (21 mm)	0.83 in (21 mm) offset 0.16 in (4 mm) from center



Excelsior OEM Controller and CDRH Power Supply Dimensions